


Protein Labeling with Pyrylium Dyes

Pyrylium dyes are fluorogenic dyes that are reactive with primary amines and form fluorescent products after conjugation with primary amino groups of peptides, proteins, and other biomolecules.



In a free state, pyrylium dyes display a weak fluorescence, but their quantum yield increases significantly after conjugation with primary amines. Also, a short-wave spectral shift of the dye fluorescence is observed. The shift in excitation/emission spectra and the increase in quantum yield after conjugation contribute to a significant reduction in the background from the unbound dye. Unbound dye is also hydrolyzed during the labeling. All this allows labeling amine-containing molecules using a simple one-step incubation at room temperature without additional purification steps.

Pyrylium dye-labeled peptides and proteins are ready to use immediately after conjugation. Therefore, they can be used in a number of «no-wash» applications, such as SDS-protein gel electrophoresis, capillary electrophoresis, isoelectric focusing, etc. Pyrylium dyes are also used as a fluorescent label in receptor binding studies.

A typical protocol for labeling proteins with pyrylium dyes is provided below.

Lumiprobe Corporation

201 International Circle, Suite 135
 Hunt Valley, Maryland 21030
 USA
 Phone: +1 888 973 6353
 Fax: +1 888 973 6354
 Email: order@lumiprobe.com

Lumiprobe GmbH

Feodor-Lynen-Strasse 23
 30625 Hannover
 Germany
 Phone: +49 511 16596811
 Fax: +49 511 16596815
 Email: de@lumiprobe.com

Lumiprobe RUS Ltd

Kotsyubinsky street, 4
 121351 Moscow
 Russian Federation
 Phone: +7 800 775 3271
 Email: ru@lumiprobe.com

Solution Preparation

1. *Dye Stock Solution*
 - a. Dissolve 1 mg of dye in 100 μ L of dimethylformamide.
Important! Do not use amine-containing solutions or buffers as a solvent.
 - b. The stock solution can be stored in the dark at 4 °C for 6 months.
2. *0.1 M Bicarbonate Buffer (pH 8.3)*
 - a. Dissolve 4.2 g NaHCO_3 in 500 mL of bidistilled water.
 - b. Adjust the buffer to pH 8.3 with 1 N NaOH.
3. *22 mM phosphate buffer (pH 7.2)*
 - a. Dissolve 5.67 g $\text{Na}_2\text{HPO}_4 \times 12\text{H}_2\text{O}$ and 0.96 g $\text{NaH}_2\text{PO}_4 \times 2\text{H}_2\text{O}$ in 1 L of bidistilled water.
 - b. Adjust the buffer to pH 7.2 with 1 N HCl.

Protein labeling

1. Dissolve 2 mg protein in 0.5 mL bicarbonate buffer (pH 8.3).
2. Add 5-10 μ L of Pyrylium-1, -4, or -6 working solution dropwise to the protein solution.
Add 10-20 μ L of the Pyrylium-8 working solution dropwise to the protein solution.
3. Gently stir the reaction mixture at room temperature for the indicated time:
 - 30 min — for Pyrylium-1;
 - 1 h — for Pyrylium-4 and -6;
 - 2 h — for Pyrylium-8.*Important!* The dye is highly reactive in the pH range from 8.0 to 9.0.
4. When incubated in the basic solution, the reaction mixture will change color from blue (Pyrylium-1), violet (Pyrylium-4), or purple (Pyrylium-6) to yellow.
5. Unbound dye is hydrolyzed during the labeling procedure, so pyrylium dye-labeled peptides and proteins are ready for use immediately after conjugation.

Purification of the conjugated protein

For some applications, purification of the dye-conjugated protein may be required.

The labeled protein is purified by size exclusion chromatography using Sephadex G25 as the stationary phase and 22 mM phosphate buffer (pH 7.2) as the eluting agent. The red band indicates the labeled protein.