

ProbeMaster® Lyo dsGreen, 5×

<http://cn.lumiprobe.com/p/probemaster-lyo-dsgreen>

ProbeMaster® Lyo dsGreen is a lyophilized master mix containing all the necessary components for performing the Polymerase Chain Reaction (PCR), as well as the intercalating dye [dsGreen](#). The composition of the mix has been optimized to achieve ideal amplification processivity and specificity. To reconstitute the mix into liquid form, simply add the specified volume of water.

The ProbeMaster® Lyo dsGreen mix is suitable for both real-time PCR using the intercalating dye dsGreen and for DNA amplification followed by result detection via electrophoresis. Because the mix lacks UDG/dUTP, it can be used for routine cloning tasks and other applications that require subsequent use of the PCR product.

Master Mix Composition

- HS Taq DNA Polymerase;
- Deoxynucleoside triphosphate (dNTP) mix;
- PCR buffer (containing Mg^{2+});
- dsGreen intercalating dye;
- Cryoprotectants

Key Features

- One tube of the lyophilized mix, once reconstituted with 450 μ L of water, is sufficient for performing 100 reactions with a volume of 25 μ L each.
- The mix is completely ready for use, which reduces the risk of sample contamination and significantly saves time during reaction setup. To set up a reaction, one needs only to add the DNA template, primers, and water to the mix.
- Suitable for PCR amplification of fragments up to 3,000 base pairs (bp) in length, with a GC content not exceeding 70%, and for applications that do not require high-fidelity amplification.
- Genomic, viral, plasmid DNA, and other types of DNA can be used as templates.
- The reaction mix contains Taq polymerase featuring «hot start» technology. The HS Taq DNA polymerase utilized in this product consists of a complex formed between monoclonal antibodies and the enzyme. Heating the sample during the initial PCR cycle inactivates the antibodies within the complex, thereby activating the enzyme. This "Hot-Start" technology effectively prevents nonspecific amplification and primer dimer formation.
- The included HS Taq DNA polymerase exhibits 5'→3' polymerase, 5'→3' exonuclease, and adenyltransferase activities, making the resulting PCR products suitable for TA cloning.
- The mixture contains dsGreen, an intercalating dye highly sensitive to the presence of double-stranded DNA, which enables «real-time» PCR (quantitative PCR) without the need for added fluorescent probes.
- Contains no UDG or dUTP.

Applications

Real-time PCR, PCR with electrophoretic detection, PCR using cDNA templates generated via reverse transcription, genotyping, colony PCR, generation of products for TA cloning, and others.

Equipment Compatibility

Compatible with all types of thermal cyclers.

PCR reaction mixture selection table

Name	Reaction mixtures for quantitative PCR (RT-PCR)				Application
	dsGreen	Eva488	ROX	UDG, dUTP	
ProbeMaster® Lyo UDG Cat.# •0514	—	—	—	✓	qPCR with DNA probes or intercalating dye
ProbeMaster® Lyo ROX Cat.# •0114	—	—	✓	—	
ProbeMaster® Lyo Eva488 Cat.# •0614	—	✓	—	—	
ProbeMaster® Lyo Eva488 ROXCat.# •0714	—	✓	✓	—	
ProbeMaster® Lyo dsGreen Cat.# •0814	✓	—	—	—	
Reaction mixture for standard PCR					
ProbeMaster® Lyo GEL Cat.# •0024	—	—	—	—	PCR followed by gel electrophoresis analysis, contains dye for application to gel
ProbeMaster® Lyo GEL UDGCat.# •0524	—	—	—	✓	
Universal reaction mixture					
ProbeMaster® Lyo UNI Cat.# •0534	—	—	—	—	qPCR with DNA probes/intercalating dye or standard PCR followed by gel electrophoresis analysis

溶
解
度

质 功能测试
量
控
制

储 储存: 4°C 下 12 个月 (自交付之日起)。

存 运输: 25°C 以下条件下 21 天。

条 复溶成液体后, 可在 4°C 下储存最多 30 天, 或在有效期内冷冻并储存在 -20°C 下。重构的混合物最多可经历五次冻融循环。

件

法 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试, 且未经明示或暗示授权用于其他任何用途, 包括但不限于体外诊断、人类或动物用途, 以及商业用途。

声
明

ProbeMaster® 是 Lumiprobe 的注册商标。