

ProbeMaster® Lyo GEL UDG, 5x

<http://cn.lumiprobe.com/p/probemaster-lyo-gel-udg>

ProbeMaster® Lyo GEL UDG is a lyophilized, ready-to-use reaction mixture containing all the necessary components for PCR with subsequent detection by electrophoresis. To reconstitute the mixture into a liquid form, add the specified amount of water. The mixture composition is optimized to obtain ideal results in terms of processivity and amplification specificity. Thanks to the high mixture density and the presence of dyes (Bromophenol Blue and Xylene Cyanol), the sample does not need to be mixed with loading buffer before application to the gel. The presence of two dyes allows for precise control of electrophoresis time.

The ProbeMaster® Lyo GEL UDG reaction mixture is suitable for DNA amplification, with results subsequently detected by electrophoresis. Uracil-DNA glycosylase eliminates contamination from amplicons generated in previous reactions and prevents false-positive results, especially during electrophoresis. Due to the presence of dUTP, the mixture is not suitable for applications where the amplification products need to be used further. For such applications, we recommend using our [ProbeMaster®Lyo GEL](#) reaction mixture.

Reaction mixture composition

- HS Taq DNA polymerase;
- Uracil-DNA glycosylase (UDG);
- Deoxynucleoside triphosphate mixture (including dUTP);
- PCR buffer (contains Mg²⁺);
- Dyes for gel loading;
- Cryoprotectants

Key characteristics

- One tube of lyophilized mixture, after dilution in 450 µL of water, is sufficient for 100 reactions of 25 µL each.
- The mixture is completely ready for use. To set up the reaction, only the DNA sample, primers, and water need to be added to the mixture, which significantly saves time. The ready-to-use reaction mixture format reduces the risk of sample contamination.
- Uracil-DNA glycosylase removes amplicon contamination from previous reactions and prevents false-positive results, which is especially important when performing amplicon electrophoresis.
- Suitable for PCR of fragments up to 3000 bp in length, with no more than 70% GC content, and not requiring high-precision amplification.
- Genomic, viral, plasmid DNA, and other templates, as well as cDNA obtained by reverse transcription, can be used.
- Contains a highly processive Hot-Start Taq polymerase, activated at 95 °C for 1 minute. The HS Taq DNA polymerase is a complex of monoclonal antibodies with the enzyme. Heating the sample in the first PCR cycle inactivates the antibodies in the complex and activates the enzyme. The "Hot-Start" technology prevents non-specific amplification and primer dimer formation.
- The composition and density of the mixture are optimized for direct application of the sample to an agarose gel after amplification.
- Due to the dyes included in the mixture, samples are easy to load onto an agarose gel. The presence of two dyes (Bromophenol Blue and Xylene Cyanol) allows for precise control of electrophoresis time.

Applications

PCR with detection of amplification products using gel electrophoresis, and PCR after reverse transcription.

Equipment compatibility

Compatible with any thermal cycler.

PCR reaction mixture selection table

Name	Reaction mixtures for quantitative PCR (RT-PCR)				Application
	dsGreen	Eva488	ROX	UDG, dUTP	
ProbeMaster[®] Lyo UDG Cat.# •0514	—	—	—	✓	
ProbeMaster[®] Lyo ROX Cat.# •0114	—	—	✓	—	
ProbeMaster[®] Lyo Eva488 Cat.# •0614	—	✓	—	—	qPCR with DNA probes or intercalating dye
ProbeMaster[®] Lyo Eva488 ROXCat.# •0714	—	✓	✓	—	
ProbeMaster[®] Lyo dsGreen Cat.# •0814	✓	—	—	—	
Reaction mixture for standard PCR					
ProbeMaster[®] Lyo GEL Cat.# •0024	—	—	—	—	PCR followed by gel electrophoresis analysis, contains dye for application to gel
ProbeMaster[®] Lyo GEL UDGCat.# •0524	—	—	—	✓	
Universal reaction mixture					
ProbeMaster[®] Lyo UNI Cat.# •0534	—	—	—	—	qPCR with DNA probes/intercalating dye or standard PCR followed by gel electrophoresis analysis

外观

溶解度

质量控制

储存: 4°C 下 12 个月 (自交付之日起)。

运输: 在 25°C 温度下最多 21 天。

重新配制为液体形式后, 可在 4°C 下储存最多 30 天, 或在有效期内冷冻并储存在 -20°C 下。重构的混合物最多可经历五次冻融循环。

件

本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试, 且未经明示或暗示授权用于其他任何用途, 包括但不限于体外诊断、人类或动物用途, 以及商业用途。

声明

ProbeMaster[®] 是 Lumiprobe 的注册商标。