

---

## LumiMount® Fluorescence Mounting Medium

<http://cn.lumiprobe.com/p/lumimount-mounting-medium>

LumiMount® is an aqueous-based mounting medium for coverslipping cells and tissues stained with various fluorescent dyes.

LumiMount is compatible with most fluorescent markers used in microscopy applications and can reduce fluorochrome quenching during imaging due to anti-fade agent content. The refractive index of liquid mounting medium is 1.38, but it can increase to 1.46 as the water evaporates at the edges of the coverslip.

LumiMount provides a semi-permanent seal for prolonged slide storage at 2°C to 8°C. Since LumiMount is water-soluble, the coverslip may be removed by submerging the slide in a PBS solution until the coverslip is loosened. For better sealing of slides, it is recommended to seal the edges of the coverslip with waterproof nail polish.

---

外观  
溶解度  
质量控制

储存条件  
储存: 收到后室温避光保存24个月。不要冷藏。避免暴露在光线下。

法律声明  
本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试, 且未经明示或暗示授权用于其他任何用途, 包括但不限于体外诊断、人类或动物用途, 以及商业用途。