
LumiMount® Plus Fluorescence Mounting Medium

<http://cn.lumiprobe.com/p/lumimount-plus-mounting-medium>

LumiMount® Plus is a ready-to-use aqueous-based mounting medium for coverslipping and imaging fluorescently labeled cell or tissue samples.

LumiMount Plus has a viscosity that provides easy application and uniform spreading over the tissue section, and it is not prone to bubbles.

Due to its unique formulation, LumiMount Plus reduces the rapid photobleaching of fluorescent proteins and dyes and prolongs possible imaging time. The refractive index of liquid LumiMount Plus is 1.43, but this property increases up to 1.52 after water evaporation and solidification of mounting.

LumiMount Plus provides a semi-permanent seal for prolonged slide storage at 2° to 8°C. Since LumiMount Plus is water-soluble, the coverslip may be removed by submerging the slide in a PBS solution until the coverslip is loosened.

LumiMount Plus is not compatible with intercalating dyes such as Hoechst, DAPI etc.

外
观
溶
解
度
质
量
控
制

储存: 收到后室温避光保存24个月。不要冷藏。避免暴露在光线下。

存
条
件

法 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试,且未经明示或暗示授权用于其他任何用途,包括但不限于体外诊
律 断、人类或动物用途,以及商业用途。
声
明