

---

## Proteinase K Dilution Buffer

<http://cn.lumiprobe.com/p/proteinase-k-dilution-buffer>

Buffer solution for dissolving lyophilized proteinase K, maintaining pH 7.5. The buffer contains calcium chloride, which increases the thermal stability of proteinase K.

The recommended concentration of proteinase K stock solution is 20 mg/mL. After dilution, the proteinase K solution must be stored at a temperature of  $-25$  to  $-18$  °C.

This buffer is included in kits for DNA isolation [Spin Kit for Blood and Buccal Swab](#) and [Spin Kit for Any Sample](#).

---

外观  
质量  
控制  
储存  
条件

本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。

声明