

## **QuDye® dsDNA HS Reagent, 200×**

<http://cn.lumiprobe.com/p/qureader-dsdna-hs-reagent>

QuDye® dsDNA HS Reagent 200×, quantification reagent is extremely selective for dsDNA over ssDNA and RNA. To prepare the working solution the concentrated QuDye dsDNA HS Reagent has to be diluted 1:200.

QuDye dsDNA HS Reagent usage is a high sensitive method over other traditional of DNA quantitation and is intended for use with all models fluorometers.

The 200× formulation is also a component of the nucleic acid quantification [QuDye dsDNA HS Assay Kit](#) for dsDNA samples measurements.

---

法 本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。  
声  
明