



## FAM phosphoramidite, 5-isomer

<http://cn.lumiprobe.com/p/5-fam-phosphoramidite>

Standard fluorescein (FAM) phosphoramidite for 5'-terminal oligonucleotide labeling, high isomeric purity single isomer. This highly purified amidite reagent ensures excellent coupling results with various synthesizers.

Samples and special offers are available for oligo manufacturers. Contact us for bulk supplies of this reagent.



外观:

分子 843.94

量:

CAS 147566-42-5

编号:

分子  $C_{46}H_{56}N_3O_{10}P$

式:

溶解

度:

质量

控制:

储存

条件:

法律 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。

激发/ 492

吸收

极大

值,

纳米:

$\epsilon$ , 摩 72000

尔吸

光系

数  $m^2$

发射 517

极大

值,

纳米:

荧光 0.93

量子

产率:

$CF_{260}$ : 0.22

$CF_{280}$ : 0.17

稀释

剂:

偶联

条件:

切割

条件:

解保

护条

件: