


Aldosterone-d4

<http://cn.lumiprobe.com/p/aldosterone-d4>

Aldosterone is a steroid hormone synthesized in the adrenal cortex. It plays a key role in regulating water-salt balance and blood pressure through its involvement in the renin-angiotensin-aldosterone system (RAAS). Aldosterone is used as a clinical marker in the diagnosis of primary aldosteronism (PA), as well as for assessing RAAS activity in renal disease, heart failure, and other conditions associated with water-electrolyte imbalance.

The product is intended for research purposes only.



外观: 白色粉末

分子量: 364.47

编号: CAS 2368815-22-7

分子式: $C_{21}H_{24}D_4O_5$

溶解度:

质量控制: NMR 1H 和HPLC-MS (95+ %, d: 98+ %)

储存条件:

收到后 -20°C 避光保存 24 个月。 运输: 室温下最多可保存3周。 干燥。

法律声明:

本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。