


## Alkyne CPG 1000A

<http://cn.lumiprobe.com/p/alkyne-cpg-1000>

Controlled pore glass solid support for synthesizing 3'-alkyne modified oligonucleotides. The support allows the preparation of oligos of up to 100mers. The modification is stable in standard deblocking conditions and can tolerate ammonia, AMA, and other basic reagents.



外观:

法律声明: 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试, 且未经明示或暗示授权用于其他任何用途, 包括但不限于体外诊断、人类或动物用途, 以及商业用途。

孔径大小, 埃: 1000

典型载 40-70  
荷  $\mu\text{mol/g}$

偶联条件:

切割条件: