


## Alkyne-PEG3-amine

<http://cn.lumiprobe.com/p/alkyne-peg3-amine-932741-19-0>

Alkyne-PEG3-amine is a heterobifunctional molecule based on triethyleneglycol (PEG3) scaffold. It contains amine and alkyne functional groups. The former can be conjugated with various electrophiles. The latter is suitable for copper-catalyzed click chemistry and Sonogashira coupling.

PEG linkers are hydrophilic linkers of choice for conjugating biomolecules.



外观:

分子 187.24

子  
量:

CAS 932741-19-0

编号:

分子 C<sub>9</sub>H<sub>17</sub>NO<sub>3</sub>

式:

溶解

度:

质量

控制:

储存

条件:

法律

声明:

本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试, 且未经明示或暗示授权用于其他任何用途, 包括但不限于体外诊断、人类或动物用途, 以及商业用途。