


Alkyne-PEG3-iodide

<http://cn.lumiprobe.com/p/alkyne-peg3-iodide>

Alkyne-PEG3-iodide (or Iodo-PEG3-alkyne) is a hydrophilic bifunctional PEG derivative for the alkylation of nucleophilic groups and copper-catalyzed click chemistry. As compared with chloride and bromide, iodide is the halogen group that is the most active in alkylation reactions.



外观:

分子量: 298.12

量:

CAS 1234387-35-9

编号:

分子式: C₉H₁₅I O₃

式:

IUPAC 3-[2-[2-(2-iodoethoxy)ethoxy]ethoxy]-1-propyne

名称:

溶解

度:

质量

控制:

储存

条件:

法律声明: 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试, 且未经明示或暗示授权用于其他任何用途, 包括但不限于体外诊断、人类或动物用途, 以及商业用途。