


Alkyne-PEG4-amine

<http://cn.lumiprobe.com/p/alkyne-peg4-amine-1013921-36-2>

Alkyne-PEG4-amine is a bifunctional oligoethylene glycol linker containing both an amine group for the easy attachment to various electrophilic reagents (NHS esters, epoxides, etc.), and an alkyne group for click chemistry. This allows a smooth conjugation of different biomolecules.

Reasonably long, highly hydrophilic PEG4 linker ensures smooth conjugation and separation of conjugated fragments in aqueous solutions.



外观:

分子 231.29

量:

CAS 1013921-36-2

编号:

分子 C₁₁H₂₁NO₄

式:

IUPAC 3,6,9,12-Tetraoxapentadec-14-yn-1-amine

名称:

溶解

度:

质量

控制:

储存

条件:

法律 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊
声明：断、人类或动物用途，以及商业用途。