


Androstenedione-d7

<http://cn.lumiprobe.com/p/androstenedione-d7>

Androstenedione is a steroid hormone produced in testes in the ovary and in the adrenal cortex. Androstenedione is a precursor of testosterone and estrone. Androstenedione originates either from the conversion of dehydroepiandrosterone or from 17-hydroxyprogesterone. Androstenedione (ASD) is a biomarker used in the diagnostic workup of hyperandrogenism, congenital adrenal hyperplasia, premature adrenarche, and polycystic ovary syndrome (PCOS).



外观:

分子 293.45

量:

CAS 67034-85-9

编号:

分子 $C_{19}H_{19}D_7O_2$

式:

质量 NMR 1H and HPLC-MS (95+ %, D: 98+ %)

控制:

储

收到后在 -20 °C 黑暗条件下可保存 24 个月。运输: 室温下最多可保存 3 周。干燥。

存

条

件:

法律 本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊

断、人类或动物用途，以及商业用途。

声明: