

---

## Annexin V, His Tag

<http://cn.lumiprobe.com/p/annexin-v>

Annexin V (or Annexin A5) belongs to the phospholipid-binding Annexin family of intracellular proteins. In flow cytometry and fluorescent microscopy, Annexin V is commonly used to detect apoptotic cells by its ability to specifically bind to phosphatidylserine, which migrates from the inner to the outer leaflet of the plasma membrane during the early stages of apoptosis.

Annexin V protein can be used for inhibiting the procoagulant and proinflammatory activities of the apoptotic cells. It also can be labeled with fluorescent dye using one of our [Protein Labeling Kits](#). This protein carries a polyhistidine tag at the C-terminus and can be recognized with Anti-His Tag Antibodies.

---

分 37.0 kDa

子  
量

质  
量  
控  
制

储  
存  
条  
件

法 本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊  
律 断、人类或动物用途，以及商业用途。

声  
明