


Azide-PEG3-amine

<http://cn.lumiprobe.com/p/azido-peg3-amine-cas-166388-57-4>

Azide-PEG3-amine is a bifunctional reagent with an azide group, and an amine group flanking triethylene glycol. This linker is useful for the conjugation of carboxylic derivatives (by acylation of the amine group) with alkynes (via a click chemistry reaction).



外观:

分子量: 174.2

CAS 编号: 166388-57-4

分子式: C₆H₁₄N₄O₂

IUPAC 名称: 1-Amino-8-azido-3,6-dioxaoctane

溶解度:

质量控制:

储存条件: