


## Biotin-11-ddUTP

<http://cn.lumiprobe.com/p/biotin-11-ddutp>

Biotin-11-ddUTP is a biotinylated dideoxyuridine triphosphate analog. Its dideoxy structure acts as an irreversible chain terminator during DNA synthesis, while the biotin moiety provides a stable tag for detection. The C11 spacer arm optimizes the molecule by facilitating polymerase recognition and ensuring unhindered interaction between the conjugated biotin and streptavidin.

Biotin-11-ddUTP can be incorporated by many DNA polymerases and terminal transferases into growing DNA strands, where it irreversibly terminates elongation while introducing a single biotin tag at the incorporation site.

This reagent is used for the enzymatic preparation of biotin-labeled oligonucleotide probes, end-labeling of DNA fragments, DNA sequencing, SNP detection, and probe labeling.



外观

分子 862.47

分子量

分子式 C<sub>28</sub>H<sub>40</sub>Li<sub>3</sub>N<sub>6</sub>O<sub>16</sub>P<sub>3</sub>S

分子式

质量控制

控制

储存

条件

件

法

本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。

律

声明

明