



Cyanine5 Boc-hydrazide

<http://cn.lumiprobe.com/p/cy5-boc-hydrazide>

This product will be discontinued soon. We may still have some material in stock (please see availability above).

Boc-protected Cyanine5 hydrazide reactive dye for the labeling of aldehydes and ketones.

Boc- (tert-butoxycarbonyl) protective group can be removed with trifluoroacetic acid, or other strong acid. Reactive hydrazide is formed, which can then couple with various carbonyl groups yielding Cyanine5-labeled molecules.



外观:

分子量: 633.26

分子式: $C_{37}H_{49}ClN_4O_3$

溶解度:

质量控制:

储存条件:

激发/吸收极大值, 纳米: 646

ϵ , 摩尔吸光系数, cm^{-1} : 250000

发射极大值, 纳米: 662

荧光量子产率: 0.2

CF_{260} : 0.03

CF_{280} : 0.04