


## Cyanine7.5 dimethyl

<http://cn.lumiprobe.com/p/cyanine75-dimethyl>

Non-functionalized Cyanine7.5 dye derivative with an aliphatic side chain and good solubility in organic solvents. The dye can be used as a non-reactive fluorophore for control experiments and for calibration.



外观:

分子量: 585.24

分子式: C<sub>40</sub>H<sub>41</sub>ClN<sub>2</sub>

溶解度:

质量控制:

储存条件:

激发/吸收极大值, 纳米: 788

ε, 摩尔吸光系数, cm<sup>-1</sup>: 223000

发射极大值, 纳米: 808

荧光量子产率: 0.10