


## Cyanine7.5 dimethyl

<http://cn.lumiprobe.com/p/cyanine75-dimethyl>

Non-functionalized Cyanine7.5 dye derivative with an aliphatic side chain and good solubility in organic solvents. The dye can be used as a non-reactive fluorophore for control experiments and for calibration.



外观:

分子量: 585.24

分子式:  $C_{40}H_{41}ClN_2$

溶解度:

质量控制:

储存条件:

激发/吸收极大值, 纳米: 788

$\epsilon$ , 摩尔吸光系数,  $cm^{-1}$ : 223000

发射极大值, 纳米: 808

荧光量子产率: 0.10