

D-Luciferin, Potassium Salt


<http://cn.lumiprobe.com/p/d-luciferin>

D-luciferin is a low-molecular-weight substrate for the firefly (*Photinus pyralis*) luciferase enzyme. The ATP-dependent oxidation of D-luciferin by luciferase is accompanied by intense chemiluminescence at 560 nm.

Luciferase is encoded by the *luc* gene, which is widely used as a reporter gene for the genetic modification of cells and organisms. When D-luciferin is administered, it is taken up by cells, and the appearance of luminescence indicates the presence and activity of the luciferase. Because most organisms do not naturally synthesize luciferase, the background signal is virtually absent, making this reporter system extremely sensitive.

The luciferase/luciferin system is used to study gene expression, noninvasively monitor tumor growth, metastasis, or the fate of transplanted stem cells in living animals over time, visualize infections, measure ATP levels in cells, and assess cell viability.

For *in vivo* experiments, the standard working dose is 100-150 mg/kg of mouse body weight administered intraperitoneally. For *in vitro* assays, the working concentration is typically 0.1-1.0 mM. The signal intensity and kinetics depend on pH, ATP and Mg^{2+} concentrations, temperature, and the permeability of cell membranes to the substrate.



外观:

分子量: 318.42

分子式:

CAS 115144-35-9 (potassium salt)

编号:

分子式:

$C_{11}H_7KN_2O_3S_2$

溶解度:

质量

控制:

储存

条件:

法律

声明:

本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。