


## dsGreen, solid form

<http://cn.lumiprobe.com/p/dsgreen-i-solid-form>

dsGreen is a highly sensitive and selective stain for double-stranded DNA (dsDNA). dsGreen solutions are used for DNA binding and [gel staining](#). This product is a solid form of dsGreen that can be used to prepare solutions according to research protocols. Lumiprobe also manufactures dsGreen in the form of [a solution](#) specifically optimized for real-time PCR.



外观:  
分子 509.73  
量:  
CAS 163795-75-3  
编号:  
分子  $C_{32}H_{37}N_4S$   
式:  
溶解  
度:  
质量  
控制:  
储存  
条件:

法律声明: 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试,且未经明示或暗示授权用于其他任何用途,包括但不限于体外诊断、人类或动物用途,以及商业用途。

激发/吸收  
极大值,  
纳米:  
发射  
极大值,  
纳米: