


dsGreen®, solid form

<http://cn.lumiprobe.com/p/dsgreen-i-solid-form>

dsGreen® is a highly sensitive and selective stain for double-stranded DNA (dsDNA). dsGreen solutions are used for DNA binding and [gel staining](#). This product is a solid form of dsGreen that can be used to prepare solutions according to research protocols. Lumiprobe also manufactures dsGreen in the form of [a solution](#) specifically optimized for real-time PCR.



外观:
分子 509.73
量:
CAS 163795-75-3
编号:
分子 C₃₂H₃₇N₄S
式:
溶解
度:
质量
控制:
储存
条件:

法律声明: 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试, 且未经明示或暗示授权用于其他任何用途, 包括但不限于体外诊断、人类或动物用途, 以及商业用途。

激发/吸收
极大值,
纳米:
发射
极大值,
纳米: