


dsGreen, solid form

<http://cn.lumiprobe.com/p/dsgreen-i-solid-form>

dsGreen is a highly sensitive and selective stain for double-stranded DNA (dsDNA). dsGreen solutions are used for DNA binding and [gel staining](#). This product is a solid form of dsGreen that can be used to prepare solutions according to research protocols. Lumiprobe also manufactures dsGreen in the form of [a solution](#) specifically optimized for real-time PCR.



外

观:

分子 509.73

量:

CAS 163795-75-3

编

号:

分子 $C_{32}H_{37}N_4S$

式:

溶解

度:

质量

控

制:

储存

条

件:

法律 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试,且未经明示或暗示授权用于其他任何用途,包括但不限于体外诊

声 明:

激 454

发/

吸收

极大

值,

纳

米:

发射 524

极大

值,

纳

米: