



Lumiprobe Corporation

201 International Circle, 135号套房

马里兰州亨特瓦雷, 21030

美国

电话: +1 888 973 6353

传真: +1 888 973 6354

电子邮件: order@lumiprobe.com

QuReader® 1 Green

<http://cn.lumiprobe.com/p/fluorometer-qureader-1-well-green>

The QuReader® 1 Green fluorometer contains one fluorescence channel and enables dsDNA, ssDNA, RNA, and protein quantitation using fluorescence measuring in green channel. Total compatibility with the QuDye® assays ensure highly sensitive determination of nucleic acids and protein concentrations in the lower pg/uL range.

法 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。
声
明