



Lumiprobe Corporation
201 International Circle, 135号套房
马里兰州亨特瓦雷, 21030
美国
电话: +1 888 973 6353
传真: +1 888 973 6354
电子邮件: order@lumiprobe.com

QuReader® 8 Green

<http://cn.lumiprobe.com/p/fluorometer-qureader-8-well-green>

The QuReader® 8 Green fluorometer is a desk fluorometer with a touch screen designed for dsDNA, ssDNA, RNA, and protein quantitation. The QuReader® fluorometer enables the fluorescence measuring for up to eight samples using green channel. Total compatibility with the QuDye® assays ensure highly sensitive determination of nucleic acid and protein concentrations in the lower pg/uL range.

法 本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。

声 明