

---


## Iodo-PEG3-carboxylic acid

<http://cn.lumiprobe.com/p/i-peg3-carboxylic-acid>

A molecule containing carboxylic acid (COOH) group and alkyl iodide functional group.

Alkyl iodides are powerful alkylators. They react with a variety of nucleophiles, very active against thiol groups, amines, alcoholates, and phenolates.

This molecule allows performing a variety of conjugations via a flexible and hydrophilic PEG3 linker.



---

外观:

分子量: 332.13

分子式: C<sub>9</sub>H<sub>17</sub>IO<sub>5</sub>

溶解度:

质量控制:

储存条件: