

Iodo-PEG3-iodide

<http://cn.lumiprobe.com/p/iodo-peg3-iodide>

PEG3 diiodide is a homobifunctional linker with two iodide functional groups which are very reactive against many nucleophiles like thiols, amines, alkoxy anions, etc. Can be used to synthesize hydrophilic symmetrical conjugates.

外观:

分子量: 369.97

CAS 编号: 36839-55-1

分子式: $C_6H_{12}I_2O_2$

IUPAC 名称: 1,2-bis(2-iodoethoxy)ethane

溶解度:

质量控制:

储存条件: