
LumiTracker® Myelin Cyan

<http://cn.lumiprobe.com/p/lumitracker-myelin-cyan>

LumiTracker® Myelin Cyan is a fluorescent myelin probe that selectively labels intact myelin sheaths in the central nervous system. The dye produces staining patterns closely matching those obtained with myelin basic protein (MBP) immunostaining and preferentially labels compact myelin in white matter tracts.

Unlike conventional myelin stains, LumiTracker® Myelin Cyan readily crosses the blood-brain barrier and accumulates in the brain following intravenous administration. This property enables direct visualization of myelinated structures *in vivo* and *ex vivo*, making the probe suitable for studies of myelin distribution and integrity.

LumiTracker® Myelin Cyan fluorescence is markedly reduced in demyelinated tissue and in myelin-deficient animal models, allowing reliable detection of demyelinating lesions. These characteristics make LumiTracker® Myelin Cyan a valuable tool for investigating myelin pathology and monitoring demyelination in experimental models of neurological disease.

外观
质量
控制
储存
条件

法 本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。

声明

LumiTracker® 是 Lumiprobe 的商标