

## **Magnetic Beads (100 mg/mL)**

<http://cn.lumiprobe.com/p/magnetic-beads>

Magnetic silica-based particles represent the modern standard for rapid and efficient DNA and RNA isolation. This technology delivers preparations of the highest purity, free from inhibitors and with minimal sample loss. Offering seamless compatibility with automated platforms, it is equally suited for manual protocols and high-throughput research. You will obtain pure, ready-to-use nucleic acids, ideal for even the most demanding downstream applications.

---

法 本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。

声  
明