

Azide-PEG3-azide

<http://cn.lumiprobe.com/p/peg3-diazide-59559-06-7>

PEG3 (triethylene glycol) diazide (1,8-diazido-3,6-dioxaoctane) is a symmetrical bifunctional crosslinker bearing two azide groups.



外观:

分子量: 200.2

CAS 编号: 59559-06-7

分子式: $C_6H_{12}N_6O_2$

IUPAC 名称: 1,8-diazido-3,6-dioxaoctane, 1,2-bis(2-azidoethoxy)ethane

溶解度:

质量控制:

储存条件: