


PEG3-dicarboxylic acid

<http://cn.lumiprobe.com/p/peg3-dicarboxylic-acid-31127-85-2>

PEG3-dicarboxylic acid is a homobifunctional linker with two carboxylic acid groups flanking the PEG3 molecule.

The linker is highly hydrophilic and useful for conjugating and immobilizing proteins, peptides, and other biomolecules. The carboxy groups can be activated by carbodiimides.



分子 294.3

量:

CAS 31127-85-2

编号:

分子 C₁₂H₂₂O₈

式:

IUPAC 3-[2-[2-(2-(2-carboxyethoxy)ethoxy)ethoxy]ethoxy]propanoic acid, 4,7,10,13-tetraoxahexadecanedioic acid

名称:

溶解

度:

质量

控制:

储存

条件:

法律 本产品仅供研究目的提供和销售。 本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。