


Phosphate CPG 1000 solid support

<http://cn.lumiprobe.com/p/phosphate-cpg-1000>

Phosphate modified oligos are used for gene synthesis by ligation, and other applications. Phosphate CPG 1000 is a modifier solid support for the synthesis of oligonucleotides bearing phosphate at 3'-end.

CPG 1000 is suitable for the synthesis of long oligos up to 100 base pairs, the use of long oligos allows for a more efficient gene assembly. The CPG is compatible with standard synthesis and deblocking procedures, it ensures smooth cleavage and complete phosphorylation.



外观:

质量控制:

制:

储存条件:

件:

法律声明: 本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试,且未经明示或暗示授权用于其他任何用途,包括但不限于体外诊断、人类或动物用途,以及商业用途。

孔径大小,埃: 1000

典型载

荷, $\mu\text{mol/g}$: 30-50

偶联条件:

件:

解保护

条件: