



## Pyrene azide 3

<http://cn.lumiprobe.com/p/pyrene-azide-3>

**This product will be discontinued soon in favor of [Pyrene azide 1](#). We may still have some material in stock (please see availability above).**

This pyrene azide is a derivative of pyrenebutyric acid. Contains triethylene glycol linker to facilitate dissolution during labeling reaction.

Other pyrene azides ([1](#), including detailed [description](#) of fluorophore, and [2](#)) are available.



外观:

分子量: 444.53

CAS 编号: 2252168-06-0

分子式: C<sub>26</sub>H<sub>28</sub>N<sub>4</sub>O<sub>3</sub>

溶解度:

质量控制:

储存条件:

激发/吸收极大值, 纳米: 343; 326; 313; 276; 265; 242; 234

发射极大值, 纳米: 377; 397