



sulfo-Cyanine7 bis-NHS ester

<http://cn.lumiprobe.com/p/sulfo-cy7-bis-nhs-ester>

sulfo-Cyanine7 bis-NHS ester is a water-soluble, bifunctional cross-linker containing a core of near-infrared, water-soluble sulfo-Cyanine7, and two amine-reactive NHS ester groups.



外观:

分子量: 1041.23

分子式: $C_{50}H_{57}N_4KO_4S_2$

溶解度:

质量控制:

储存条件:

法律声明:

激发/吸收极大值, 纳米: 750

ϵ , 摩尔吸光系数, cm^{-1} : 240600

发射极大值, 纳米: 773

荧光量子产率: 0.24

CF_{260} : 0.04

CF_{280} : 0.04