


TETA-PEG4-azide hydrochloride

<http://cn.lumiprobe.com/p/teta-peg4-azide-hydrochloride>

TETA (1,4,8,11-Tetraazacyclotetradecane-1,4,8,11-tetraacetic acid) is a complexing agent for di- and trivalent metal ions. TETA-PEG4-azide is a bifunctional linker containing chelating and azide groups for conjugation with alkyne-containing biomolecules. The structure of the reagent features a long hydrophilic PEG4 linker that separates the TETA residue from the target molecule. The linker also enhances aqueous solubility to facilitate conjugation.



外观

分子 778.56

分子量

分子式 $C_{26}H_{52}N_8Cl_4O_{10}$

分子式

溶解度

质量控制

储存条件

法律声明
本产品仅供研究目的提供和销售。本产品并未经过食品、药品、医疗器械、化妆品等领域的安全性和效力测试，且未经明示或暗示授权用于其他任何用途，包括但不限于体外诊断、人类或动物用途，以及商业用途。